



## XTI program to view the SXR tomography data

- The XTI program, written in IDL, can be used either on dave or aurora/juno. In dave log on with username *tomography*; in aurora/juno log on as *mstdata*. Of course, in aurora/juno you will see only the shots of a couple of days.
- In dave type:
  - idl
  - xti, /mst
- In aurora/juno type:
  - fd
  - idl
  - xti, /mst
- The main program will start:



- To visualize the SXR data (brightness) for a shot type the following parameters:
  - 1060418098 in the Shot field
  - 14 in the t Start (ms) field
  - 21 in the t End (ms) field
  - 0.25 in the Delta (ms) field
- Click on the OK button. The raw data will be read, together with the magnetic data (Read magnetics button selected).


- To plot the four SXR profiles click on the ALL profiles button in the top-right of the window:



Move the time slider on the bottom to plot the profiles for different time or type a value in the time field, followed by enter, to plot the profiles for that time (or for the closest time available). In this case the two profiles SXR1 and SXR2 have the same Y scale, while SXR3 and SXR4 have the same scale but different form that of SXR1,SXR2. In fact, SXR1 and SXR2 have 303 microns of Beryllium, while SXR3,SXR4 have 761 microns.


- For this example the data are read in the range 14-21 ms, averaged by 0.25 ms (then the time base will be 14.00, 14.25, 14.50, 14.75, ... ms). In general, using smaller Delta will give signals with higher fluctuations (and noise); using Delta of 0.2 or 0.5 ms will give smoother signals and profiles, but will cancel also the  $m=1$  rotating structure. The sampling frequency is 500kHz, bandwidth of the electronics is about 100kHz. The raw data are in W/m<sup>2</sup>; amplifier gains, geometrical factors (solid angles of the lines of sight) etc. are all automatically calculated and included in the calculations.

- To plot the SXR and magnetic data (toroidal measurements) click on the SXR and magnetics



button:

You can select the time for the SXR profiles using the time slider or field but you can also select a time in the Bt (mT) plot (bottom-right) by clicking with the mouse (left button). You can also zoom



by selecting a region in the same plot:

To replot click on the **Replot** button or click in any point with the right button of the mouse. In this plot the following signals are plotted:

- left column, top to bottom: the four contour plot of the four SXR cameras (X axis=time, Y axis=radius) and the SXR profiles for the selected time (see the field **time**), on the left SXR1 (red) and SXR2 (white), on the right SXR3 (green) and SXR4 (white).
- right column, top to bottom: the  $m=1$  spectrum in the selected time, mode velocity, poloidal position of the mode at the toroidal position of the tomography and the mode amplitude as a function of time. For these three graphs the various  $n$  number are color coded (see the legend in the small window below the main plot window).